

Technical Manual

ChromaDazzle Glucose Dehydrogenase Activity Assay Kit

Catalogue Code: BA0024

Pack Size: 100 assays

Research Use Only

DESCRIPTION

GLUCOSE DEHYDROGENASE (GDH) belongs to the family of oxio-reductases, specifically those acting on the CH-OH group of donor with other acceptors. GDH participates in the pentose phosphate pathway. The Assay Genie non-radioactive, colorimetric ChromaDazzle Glucose Dehydrogenase Activity Assay Kit (GDH) assay is based on the reduction of the tetrazolium salt MTT in a NADH-coupled enzymatic reaction to a reduced form of MTT which exhibits an absorption maximum at 565 nm. The increase in absorbance at 565 nm is proportional to the enzyme activity.

KEY FEATURES

Fast and sensitive. Linear detection range (20 μ L sample): 0.5 to 200 U/L for 15 min reaction.

Convenient and high-throughput. Homogeneous "mix-incubate-measure" type assay. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

APPLICATIONS

GDH activity determination in biological samples (e.g. plasma, serum, tissue and culture media.)

KIT CONTENTS (100 TESTS IN 96-WELL PLATES)

Assay Buffer:	10 mL	Diaphorase:	120 μ L
NAD/MTT:	1 mL	Calibrator:	1.5 mL
Substrate:	1 mL		

Storage conditions. The kit is shipped at ambient temperature. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.

Precautions: reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

PROCEDURES

This assay is based on a kinetic reaction. To ensure identical incubation time, addition of Working Reagent to samples should be quick and mixing should be brief but thorough. Use of a multi-channel pipettor is recommended. Assays can be executed at any desired temperature (e.g. 25°C or 37°C).

Sample Preparation: Serum and plasma are assayed directly.

Tissue: prior to dissection, rinse tissue in phosphate buffered saline (pH 7.4) to remove blood. Homogenize tissue (50 mg) in ~ 200 μ L buffer containing 50 mM potassium phosphate (pH 7.5). Centrifuge at $10,000 \times g$ for 15 min at 4°C . Remove supernatant for assay.

Cell Lysate: collect cells by centrifugation at $2,000 \times g$ for 5 min at 4°C . For adherent cells, do not harvest cells using proteolytic enzymes; rather use a rubber policeman. Homogenize or sonicate cells in an appropriate volume of cold buffer containing 50 mM potassium phosphate (pH 7.5). Centrifuge at $10,000 \times g$ for 15 min at 4°C . Remove supernatant for assay.

All samples can be stored at -20 to -80°C for at least one month.

Reagent Preparation: Equilibrate reagents to desired reaction temperature (e.g. 25°C or 37°C). Briefly centrifuge tubes before use.

Prepare enough Working Reagent (WR) for all assay wells by mixing, for each 96-well assay: 8 μ L Substrate, 8 μ L NAD/MTT Solution, 1 μ L Diaphorase and 70 μ L Assay Buffer.

Reaction Preparation:

1. Transfer 100 μ L H_2O ($\text{OD}_{\text{H}_2\text{O}}$) and 100 μ L Calibrator (OD_{CAL}) solution into wells of a clear flat bottom 96-well plate.

- Transfer 20 μL of each sample into separate wells and then add 80 μL WR to each sample well. Tap plate briefly to mix.
- Read $\text{OD}_{565\text{nm}}$ (OD_0), and again after 15 min (OD_{15}) on a plate reader.

CALCULATION

Subtract the OD_0 from OD_{15} for each sample to compute the ΔOD_s values. GDH activity can then be calculated as follows:

$$\text{GDH Activity} = \frac{\Delta\text{OD}_s}{\varepsilon_{\text{mtt}} \cdot l} \times \frac{\text{Reaction Vol } (\mu\text{L})}{t \text{ (min)} \cdot \text{Sample Vol } (\mu\text{L})} \times n$$

$$= \frac{\Delta\text{OD}_s}{\text{OD}_{\text{CAL}} - \text{OD}_{\text{H}_2\text{O}}} \times \frac{273}{t \text{ (min)}} \times n \quad (\text{U/L})$$


where ε_{mtt} is the molar absorption coefficient of reduced MTT. l is the light pathlength which is calculated from the calibrator. OD_{CAL} and $\text{OD}_{\text{H}_2\text{O}}$ are $\text{OD}_{565\text{nm}}$ (OD_0) values of the Calibrator and water. t is the reaction time (15 min is the recommended time). Reaction Vol and Sample Vol are 100 μL and 20 μL , respectively. n is the dilution factor.

Unit definition: 1 Unit (U) of GDH will catalyze the conversion of 1 μmole of NAD to NADH per min at pH 8.2.

Note: If sample GDH activity exceeds 200 U/L, either use a shorter reaction time or dilute samples in water and repeat the assay. For samples with GDH activity < 5 U/L, the incubation time can be extended up to 2 hours.

MATERIALS REQUIRED, BUT NOT PROVIDED

Pipetting devices and accessories (e.g. multi-channel pipettor), clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), centrifuge tubes and plate reader.



LITERATURE

- Bak, TG (1967) "Studies on glucose dehydrogenase of *Aspergillus oryzae*. II Purification and physical and chemical properties". *Biochim. Biophys. Acta.* 139: 277–93.
- Brink, NG, et al. (1953) "Beef liver glucose dehydrogenase. 1. Purification and properties". *Acta Chem. Scand.* 7: 1081–1089.
- Thompson RE, Carper WR (1970) "Glucose dehydrogenase from pig liver. I. Isolation and purification". *Biochim. Biophys. Acta* 198 (3): 397–406.

Contact Details

Dublin, Ireland

Email: info@assaygenie.com

Web: www.assaygenie.com

Technical Support: Techsupport@assaygenie.com