

Technical Manual

ColorFluor Xanthine Assay Kit

Catalogue Code: BA0154

Pack Size: 100 assays

Research Use Only

3. **Working Reagent.** Prepare bulk working reagent by mixing 90 μL Assay Buffer, 1 μL XO Enzyme, 1 μL HRP Enzyme (*vortex briefly before pipetting*), and 1 μL Dye Reagent per reaction well in a clean tube. Transfer 90 μL Working Reagent into each reaction well. Tap plate to mix.
4. Incubate 30 min at room temperature, and then read optical density at 570 nm (550-585 nm) (OD_{30}).

FLUORIMETRIC PROCEDURE

For fluorimetric assays, the linear detection range is 3 to 250 μM xanthine. Dilute the standards from *Colorimetric Procedure* 10 \times with dH_2O to obtain standards at 200, 120, 60 and 0 μM Xanthine.

Transfer 10 μL standards and 10 μL samples into separate wells of a *black* 96-well plate.

Add 90 μL Working Reagent (see *Colorimetric Procedure*), tap plate to mix.

Incubate 30 min at room temperature, and then read fluorescence at $\lambda_{\text{ex/em}} = 530/585 \text{ nm}$ (F_{30}).


CALCULATION

Subtract blank OD_{30} or F_{30} (water, #4) from all standards and samples OD_{30} or F_{30} values and plot the ΔOD or ΔF against standard concentrations. Calculate the concentration using the equation below:


$$[\text{Xanthine}] = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{\text{Slope } (\mu\text{M}^{-1})} \times n \quad (\mu\text{M})$$

Where R_{Sample} and R_{Blank} are the optical density or fluorescent values of the sample and blank, respectively. *Slope* is the slope of the standard curve and n is the dilution factor.

Notes: If the calculated sample xanthine concentration is higher than 2 mM in colorimetric assay or 200 μM in fluorimetric assay, dilute sample in water and repeat the assay. Multiply result by the dilution factor (n).



96-well colorimetric assay



96-well fluorimetric standard

MATERIALS REQUIRED, BUT NOT PROVIDED

Pipetting devices, centrifuge tubes, clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), black 96-well plates (e.g. Greiner Bio-One, cat# 655900) and plate reader capable of either measuring absorbance between 550-585 nm or fluorescence intensity at $\lambda_{\text{ex/em}} = 530/585 \text{ nm}$.

LITERATURE

1. Ukena D, Schudt C, Sybrecht GW (1993). Adenosine receptor-blocking xanthines as inhibitors of phosphodiesterase isozymes. *Biochem. Pharm.* 45(4): 847–51.
2. Daly JW, Hide I, Müller CE, Shamim M (1991). Caffeine analogs: structure activity relationships at adenosine receptors. *Pharmacology.* 42(6): 309–21.
3. González MP, Terán C, Teijeira M (2008). Search for new antagonist ligands for adenosine receptors from QSAR point of view. How close are we?. *Medicinal Research Reviews.* 28(3): 329–71.

Contact Details

Dublin, Ireland

Email: info@assaygenie.com

Web: www.assaygenie.com

Technical Support: Techsupport@assaygenie.com