



Recombinant Protein Technical Manual
Recombinant Mouse SCPEP1/RISC Protein (His Tag)
RPES2245

Product Data:

Product SKU: RPES2245

Size: 20µg

Species: Mouse

Expression host: HEK293 Cells

Uniprot: NP_083299.3

Protein Information:

Molecular Mass: 49.3 kDa

AP Molecular Mass:

Tag: C-His

Bio-activity:

Purity: > 95 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: 2410018F01Rik;4833411K15Rik;AI957256;Risc

Immunogen Information:

Sequence: Met1-Glu452

Background: