

HDFP1096

Product Information

Product SKU:	HDFP1096	Expression Host:	HEK293	Size:	10µg
Target:	NPY5R	Tag:	C-Flag&Strep Tag		

Additional Information

Conjugate:	Unconjugated	Uniprot ID:	Q15761
Molecular Weight:	The human full length NPY5R-Strep protein has a MW of 50.7 kDa		

Protein Information

Background: The protein encoded by this gene is a receptor for neuropeptide Y and peptide YY. The encoded protein appears to be involved in regulating food intake, with defects in this gene being associated with eating disorders. Also, the encoded protein is involved in a pathway that protects neuroblastoma cells from chemotherapy-induced cell death, providing a possible therapeutic target against neuroblastoma. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Nov 2015]

Synonyms: NPY5-R, NPYR5, NPYY5-R

Protein Description: Human NPY5R-Strep full length protein-synthetic nanodisc

Formulation: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.

Protein Pathways: GPCRDB Class A Rhodopsin-like, Peptide GPCRs, Neurotrophin & Receptors.

Protein Families: GPCR, Transmembrane, Druggable Genome.

Usage: Research use only

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.