

HDFP929

Product Information

Product SKU:	HDFP929	Expression Host:	HEK293	Size:	10µg
Target:	AGRF4	Tag:	C-Flag&Strep Tag		

Additional Information

Conjugate:	Unconjugated	Uniprot ID:	Q8IZF3
Molecular Weight:	The human full length AGRF4-Strep protein has a MW of 77.7 kDa		

Protein Information

Background: Sequence analysis of this gene suggests that it encodes a member of the superfamily of G protein-coupled receptors. G protein-coupled receptors typically contain seven hydrophobic transmembrane domains, interact with guanine nucleotide binding regulatory proteins, and detect molecules outside the cell and act to transduce these signals into intracellular responses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2016]

Synonyms: GPR115, PGR18

Protein Description: Human AGRF4-Strep full length protein-synthetic nanodisc

Formulation: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.

Protein Pathways: -

Protein Families: Transmembrane, Druggable Genome.

Usage: Research use only

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.