

# Recombinant Human Basigin/EMMPRIN/CD147 Protein RPCB0259



## Product Information

<b>Product SKU:</b>	RPCB0259	<b>Gene ID:</b>	682	<b>Size:</b>	20µg
<b>Tag:</b>	C-hFc&His	<b>Reactivity:</b>	Human		

## Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	P35613
<b>Purity:</b>	> 95% by SDS-PAGE.		

## Protein Information

**Background:** This protein is a plasma membrane protein that is important in spermatogenesis, embryo implantation, neural network formation, and tumor progression. The encoded protein is also a member of the immunoglobulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Basigin/EMMPRIN/CD147 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.