

RPCB0283

---

## Product Information

<b>Product SKU:</b>	RPCB0283	<b>Gene ID:</b>	1947	<b>Size:</b>	20µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	P98172
<b>Purity:</b>	> 95% by SDS-PAGE.		

---

## Protein Information

**Background:** Ephrin-B1 also known as EFNB1, is a member of the ephrin family. The transmembrane-associated ephrin ligands and their Eph family of receptor tyrosine kinases are expressed by cells of the SVZ. Eph/ephrin interactions are implicated in axon guidance, neural crest cell migration, establishment of segmental boundaries, and formation of angiogenic capillary plexi. Eph receptors and ephrins are divided into two subclasses, A and B, based on binding specificities. Ephrin subclasses are further distinguished by their mode of attachment to the plasma membrane: ephrin-A ligands bind EphA receptors and are anchored to the plasma membrane via a glycosylphosphatidylinositol (GPI) linkage, whereas ephrin-B ligands bind EphB receptors and are anchored via a transmembrane domain. An exception is the EphA4 receptor, which binds both subclasses of ephrins.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Ephrin-B1/EFNB1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.

**Storage:**

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.