

**RPCB0344**

---

## Product Information

**Product SKU:** RPCB0344      **Gene ID:** 965      **Size:** 20µg  
**Tag:** C-hFc&His      **Reactivity:** Human

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** AAH05930.1  
**Purity:** > 95% by SDS-PAGE.

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human LFA-3/CD58 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.