

RPCB0394

Product Information

| | | | | | |
|---------------------|----------|--------------------|-------|--------------|------|
| Product SKU: | RPCB0394 | Gene ID: | 1118 | Size: | 20µg |
| Tag: | C-His | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------------|--------------------|-------------------|------------|
| Expression Host: | - | Swissprot: | AAI03696.1 |
| Purity: | > 95% by SDS-PAGE. | | |

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Chitinase-1/CHIT1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.01 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.