

Recombinant Human PTPRC/CD45RABC/CD45 Protein

RPCB0470



Product Information

Product SKU:	RPCB0470	Gene ID:	5788	Size:	100µg
Tag:	C-hFc	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P08575-3
Purity:	> 95% by SDS-PAGE ; > 95% by HPLC		

Protein Information

Background: PTPRC (also known as CD45), T cells require the protein tyrosine phosphatase CD45 to detect and respond to antigen because it activates the Src family kinase Lck, which phosphorylates the T cell antigen receptor (TCR) complex. CD45 activates Lck by opposing the negative regulatory kinase Csk. Paradoxically, CD45 has also been implicated in suppressing TCR signaling by dephosphorylating the same signaling motifs within the TCR complex upon which Lck acts.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human PTPRC/CD45RABC/CD45 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.