

**RPCB0682**

---

## Product Information

<b>Product SKU:</b>	RPCB0682	<b>Gene ID:</b>	100133941	<b>Size:</b>	100µg
<b>Tag:</b>	C-His&Avi	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Expression Host:</b>	-	<b>Swissprot:</b>	P25063-1
<b>Purity:</b>	> 95% by SDS-PAGE ; > 95% by HPLC		

---

## Protein Information

**Background:** CD24 is a sialoglycoprotein expressed at the surface of most B lymphocytes and differentiating neuroblasts. It is also expressed on neutrophils and neutrophil precursors from the myelocyte stage onwards. The potential for targeting CD24 in cancer therapy seems promising, as CD24 is overexpressed in many human cancers.

**Protein Description:** High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human CD24 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method

**Formulation:** Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.