

RPCB0781

Product Information

Product SKU: RPCB0781 **Gene ID:** 11006 **Size:** 10µg
Tag: C-hFc&His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q8NHJ6
Purity: > 97% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human LILRB4/ILT-3/CD85k Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.