

# Recombinant Cynomolgus Siglec-15/CD33L3 Protein

RPCB0860



## Product Information

**Product SKU:** RPCB0860      **Gene ID:** 102126346      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Cynomolgus

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** A0A2K5UY47  
**Purity:** > 95% by SDS-PAGE.

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus Siglec-15/CD33L3 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** Please contact us for more information.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, 0.1M Glycine, pH 9.0.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.