

**RPCB0861**

---

## Product Information

**Product SKU:** RPCB0861      **Gene ID:** 17395      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Mouse

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** P41245  
**Purity:** > 97% by SDS-PAGE.

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Mouse MMP-9 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of 50mM Tris, 10mM CaCl<sub>2</sub>, 150mM NaCl, pH 7.5. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.