

RPCB0875

---

## Product Information

**Product SKU:** RPCB0875      **Gene ID:** 5972      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Human

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** P00797  
**Purity:** > 97% by SDS-PAGE.

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Renin/REN Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.