

# Recombinant Human Fibrinogen like 1/FGL1 Protein

## RPCB0900



### Product Information

**Product SKU:** RPCB0900      **Gene ID:** 2267      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Human

### Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** Q08830-1  
**Purity:** > 90% by SDS-PAGE.

### Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Fibrinogen like 1/FGL1 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** Please contact us for more information.

**Formulation:** Supplied as a 0.22 µm filtered solution in 20mM Tris,300mM NaCl,10% Glycerol,pH 8.0Contact us for customized product form or formulation.

**Storage:** Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.