

Recombinant SARS-CoV-2 Spike S1+S2 ECD (S-ECD) (D614G) Protein

RPCB0924



Product Information

Product SKU: RPCB0924 **Gene ID:** 43740568 **Size:** 100µg
Tag: C-His **Reactivity:** SARS-CoV-2

Additional Information

Expression Host: HEK293 cells **Swissprot:** P0DTC2
Purity: > 95% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant SARS-CoV-2 Spike S1+S2 ECD(S-ECD)(D614G) Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Supplied as a 0.22 µm filtered solution in PBS, pH 7.4.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.