

Product Information

Product SKU:	RPCB1041	Gene ID:	12215	Size:	10µg
Tag:	C-His&Avi	Reactivity:	Mouse		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P18572-2
Purity:	> 90% by SDS-PAGE.		

Protein Information

Background: This protein is a plasma membrane protein that is important in spermatogenesis, embryo implantation, neural network formation, and tumor progression. The encoded protein is also a member of the immunoglobulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse Basigin/EMMPRIN/CD147 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.