

RPCB1085

Product Information

Product SKU:	RPCB1085	Gene ID:	15894	Size:	10µg
Tag:	C-hFc&His	Reactivity:	Mouse		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	Q3U8M7
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse ICAM-1/CD54 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.