

RPCB1088

Product Information

Product SKU: RPCB1088 **Gene ID:** 928965 **Size:** 10µg
Tag: C-His **Reactivity:** Monkeypox virus

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q8V4U9
Purity: > 95% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Monkeypox virus A30L Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.