

RPCB1134

---

## Product Information

**Product SKU:** RPCB1134      **Gene ID:** 929030      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Monkeypox virus

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** Q5QC89  
**Purity:** > 95% by SDS-PAGE.

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Monkeypox virus D6L Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** <0.1EU/µg

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.