

RPCB1200

Product Information

| | | | | | |
|---------------------|----------|--------------------|-------|--------------|------|
| Product SKU: | RPCB1200 | Gene ID: | 7031 | Size: | 10µg |
| Tag: | C-His | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------------|--------------------|-------------------|--------|
| Expression Host: | HEK293 cells | Swissprot: | P04155 |
| Purity: | > 97% by SDS-PAGE. | | |

Protein Information

Background: Trefoil Factor 1 (TFF1), also known as pS2, is one of three structurally related secreted proteins that contain trefoil domains. These domains adopt a three-leaved conformation held together by conserved intrachain disulfide bonds. TFF1 is an approximately 7 kDa peptide that plays an important role in epithelial regeneration and wound healing (1). Mature human TFF1 shares 67% amino acid sequence identity with mouse and rat TFF1. It is expressed by goblet cells of the gastric and intestinal mucosa and by conjunctival goblet cells (2-5). TFF1 is a copper-binding protein that can form disulfide-linked homodimers, associate into disulfide-linked complexes with Gastrokine 2, and form non-covalent complexes with the mucin MUC5AC (4, 6-8). Copper enhances TFF1 homodimerization as well as its ability to promote epithelial cell motility, wound healing, and the colonization of *H. pylori* in stomach and colon epithelia (9, 10). TFF1 is down-regulated during the progression from gastritis to gastric dysplasia to gastric cancer, although it is up-regulated in breast and prostate cancers (11-13). TFF1 inhibits the formation of calcium oxalate crystals, and its excretion in the urine is reduced in patients with kidney stones (14).

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human BCEI/PS2/TFF1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.001 EU/µg

Formulation:

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage:

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.