

# Recombinant Human FGFR-2 beta (IIIb)/KGFR/CD332 Protein

RPCB1598



## Product Information

**Product SKU:** RPCB1598      **Gene ID:** 2263      **Size:** 100µg  
**Tag:** C-His      **Reactivity:** Human

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** P21802-3  
**Purity:** > 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC

## Protein Information

**Background:** Four distinct genes encoding closely related FGF receptors, FGF R1 - 4, are known. All four genes for FGF Rs encode proteins with an N-terminal signal peptide, three immunoglobulin (Ig)-like domains, an acid-box region containing a run of acidic residues between the IgI and IgII domains, a transmembrane domain and the split tyrosine-kinase domain. Multiple forms of FGF R1 - 3 are generated by alternative splicing of the mRNAs. A frequent splicing event involving FGF R1 and 2 results in receptors containing all three Ig domains, referred to as the alpha? isoform, or only IgII and IgIII, referred to as the beta? isoform.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human FGFR-2 beta (IIIb)/KGFR/CD332 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** Less than 1EU per µg by the LAL method.

**Formulation:** -

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.