

RPCB1620

---

## Product Information

**Product SKU:** RPCB1620      **Gene ID:** 54583      **Size:** 100µg  
**Tag:** 6×His tag at the N-terminus      **Reactivity:** Human

---

## Additional Information

**Expression Host:** E. coli      **Swissprot:** Q9GZT9-1  
**Purity:** > 95% by SDS-PAGE.

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human EGLN1/PHD2 Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** <0.1EU/µg

**Formulation:** Supplied as a 0.22 µm filtered solution in 20mM HEPES, 250mM NaCl, pH 7.5.

**Storage:** Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.