

RPCB1973

Product Information

Product SKU:	RPCB1973	Gene ID:	20304	Size:	100µg
Tag:	N-hFc	Reactivity:	Mouse		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P30882
Purity:	> 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC		

Protein Information

Background: The CCR5 and the CCL5 ligand have been detected in some hematological malignancies, lymphomas, and a great number of solid tumors, but extensive studies on the role of the CCL5/CCR axis were performed only in a limited number of cancers.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse CCL5/RANTES Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.