

# Recombinant Human Integrin alpha V beta 5 (ITGAV&ITGB5) Protein

RPCB2068



## Product Information

**Product SKU:** RPCB2068      **Gene ID:** 3685&3693      **Size:** 100µg  
**Tag:** C-His&Avi      **Reactivity:** Human

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** P06756&P18084  
**Purity:** > 95% by Tris-Bis PAGE;> 95% by SEC-HPLC

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Integrin alpha V beta 5 (ITGAV&ITGB5) Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.