

**HDBS0177**

---

## Product Information

**Product SKU:** HDBS0177      **Host Species:** Homo sapiens      **Size:** 100 µg  
**Target:** B7-H4      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** Q7Z7D3  
**Isotype:** IgG1

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** VTCN1

**Common Name:** SGNB7H4V,SGN-B7H4V,SGNB 7H4V

**Recommended Dilutions:** Flow Cyt 1:100

**Applications:** Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.