

**HDBS0199**

---

## Product Information

**Product SKU:** HDBS0199      **Host Species:** Chimeric Humanized      **Size:** 100 µg  
**Target:** TEM1      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** Q9HCU0  
**Isotype:** IgG1(E239D,M241L,S286F)

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD248; CD164L1

**Common Name:** MORAb-004

**Recommended Dilutions:** ELISA 1:5000-10000, Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.