

HDBS0211

---

**Product Information**

**Product SKU:** HDBS0211      **Host Species:** Humanized      **Size:** 100 µg  
**Target:** NKG2D      **Reactivity:** Human

---

**Additional Information**

**Conjugate:** Unconjugated      **Uniprot:** P26718  
**Isotype:** IgG1

---

**Protein Information**

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** NKG2D;CD314;KLRK1;NK cell receptor D

**Common Name:** 5C5

**Recommended Dilutions:** ELISA 1:5000-10000, Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.