

## HDBS0315

---

### Product Information

<b>Product SKU:</b> HDBS0315	<b>Host Species:</b> Homo sapiens	<b>Size:</b> 100 µg
<b>Target:</b> IFNAR1	<b>Reactivity:</b> Human	

---

### Additional Information

<b>Conjugate:</b> Biotinylated	<b>Uniprot:</b> P17181
<b>Isotype:</b> IgG1	

---

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** IFN-R-1;CRF2-1;IFNAR

**Common Name:** MDX-1333, MEDI-546, anifrolumab-fnia

**Recommended Dilutions:** ELISA 1:5000-10000; Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.