

# Biotinylated Anti-TM4SF1 (AGX101 biosimilar) mAb

HDBS0351



## Product Information

**Product SKU:** HDBS0351      **Host Species:** Humanized      **Size:** 100 µg  
**Target:** TM4SF1      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P30408  
**Isotype:** IgG1

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** M3S1;TAAL6

**Common Name:** AGX101

**Recommended Dilutions:** Flow Cyt 1:100

**Applications:** Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.