

HDBS0430

---

**Product Information**

**Product SKU:** HDBS0430      **Host Species:** Chimeric Humanized      **Size:** 100 µg  
**Target:** IL1RL2      **Reactivity:** Human

---

**Additional Information**

**Conjugate:** Unconjugated      **Uniprot:** Q9HB29  
**Isotype:** IgG1(K97R,L117A,L118A)

---

**Protein Information**

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** ILRL2; IL-36R; IL1RRP2; IL-36R

**Common Name:** BI-655130,BI655130,BI 655130,spesolimab-sbzo

**Recommended Dilutions:** Flow Cyt 1:100

**Applications:** Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.