

**HDBS0475**

---

## Product Information

<b>Product SKU:</b>	HDBS0475	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	IL7RA	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	P16871
<b>Isotype:</b>	IgG1		

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** ILRA;CD127;IL7RA;CDW127;IMD104;sIL-7R;Inc-IL7R;IL7Ralpha;IL-7Ralpha;IL-7R-alpha

**Common Name:** ADX-914;ADX914

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.