

HDBS0531

---

## Product Information

**Product SKU:** HDBS0531      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** TMEM219      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** Q86XT9  
**Isotype:** IgG4

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** IGFBP3R;IGFBP-3R

**Common Name:** Ent-001;Ent001;Ent 001

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.