

HDBS0569

---

## Product Information

<b>Product SKU:</b>	HDBS0569	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	IL1B	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	P01584
<b>Isotype:</b>	IgG2		

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** IL-1;IL1F2;IL1beta;IL1-BETA

**Common Name:** S-78989;S78989;S 78989;VPM 087;VPM087;VPM-087;VPM87;VPM-87;XMA-005.2;XMA-0052;XMA52;XOMA-052;XOMA052;XOMA52;XOMA-52

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.