

HDBS0599

---

## Product Information

**Product SKU:** HDBS0599      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** IL6R      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** P08887  
**Isotype:** IgG1

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** IL6Q;gp80;CD126;HIES5;IL-6R;IL6RA;IL6RQ;IL-1Ra;IL-6RA;IL6QTL;IL-6R-1

**Common Name:** BCD 089;BCD-089;BCD089;BCD89;BCD-89;BCD 89;Ilsira

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.