

HDBS0706

Product Information

Product SKU: HDBS0706	Host Species: Whole mAb	Size: 100 µg
Target: CFB	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: P00751
Isotype: IgG1	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: BF;FB;BFD;GBG;CFAB;CFBD;PBF2;AHUS4;FBI12;H2-Bf;ARMD14

Common Name: NM-8074;NM 8074;NM8074

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.