

HDBS0840

Product Information

Product SKU:	HDBS0840	Host Species:	Whole mAb	Size:	100 µg
Target:	FCGRT	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P55899
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: FCRN;FcgammaRn;alpha-chain

Common Name: HBM-9161;HBM9161;HL-161;HL161;HL-61;HL61;IMVT 1401;IMVT-1401;IMVT1401;RVT 1401;RVT-1401;RVT1401

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.