

HDBS1062

Product Information

Product SKU: HDBS1062 **Host Species:** Whole mAb **Size:** 100 µg
Target: PD-1 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q15116
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

Common Name: AEX-1188;AEX1188;AEX 1188;INCMGA-0012;INCMGA 0012;INCMGA0012;INCMGA12;INCMGA-12;INCMGA 12;MGA-012;MGA012;MGA 012;MGA-12;MGA12;MGA 12;Zynyz

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.