

HDBS1078

Product Information

Product SKU: HDBS1078 **Host Species:** Whole mAb **Size:** 100 µg
Target: IL6R **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P08887
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL6Q;gp80;CD126;HIES5;IL-6R;IL6RA;IL6RQ;IL-1Ra;IL-6RA;IL6QTL;IL-6R-1

Common Name: Kevzara;REGN-88;REGN88;REGN 88;SAR-153191;SAR153191;SAR 153191

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.