

Biotinylated Anti-SOST (Setrusumab biosimilar) mAb

HDBS1082



Product Information

Product SKU: HDBS1082 **Host Species:** Whole mAb **Size:** 100 µg
Target: SOST **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q9BQB4
Isotype: IgG2

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CDD;VBCH;DAND6;SOST1

Common Name: BPS-804;BPS804;BPS 804;MOR-05813;MOR05813;MOR 05813;MOR-5813;MOR5813;MOR 5813

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.