

HDBS1135

---

## Product Information

**Product SKU:** HDBS1135      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** IL6      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P05231  
**Isotype:** IgG2

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CDF;HGF;HSF;BSF2;IL-6;BSF-2;IFNB2;IFN-beta-2

**Common Name:** RG-6179;RG6179;RG 6179;RO-7200220;RO 7200220;RO7200220

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.