

HDBS1146

Product Information

Product SKU: HDBS1146 **Host Species:** Whole mAb **Size:** 100 µg
Target: B7-H3 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q5ZPR3
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: B7H3;B7-H3;B7RP-2;4Ig-B7-H3;CD276

Common Name: MGA017;MGA-017;MGA17;MGA-17;MGA 017;MGA 17;MGA018;MGA-018;MGA18;MGA-18;MGA 018;MGA 18;AEX4089;AEX-4089;AEX 4089

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.