

**GenieFluor Red 780 Anti-Rat CD4
(domain 1) Antibody [OX-38]**
Catalogue Code: AGEL3117

Antibody Data

Product SKU:	AGEL3117	Clone:	OX-38
Applications:	FCM		
Reactivity:	Rat		


Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names: T-cell surface glycoprotein CD4;CD4;T-cell surface antigen T4/Leu-3;CD4;
Uniprot ID: P05540
Background: CD4, also known as T4, is a 55kD glycoprotein member of the immunoglobulin superfamily and is expressed on majority of thymocytes, macrophages, and a peripheral T cell subset (T helper cells). CD4 is a T cell co-receptor that interacts with the MHC class II molecule and is involved in T cell activation.

Form: Liquid
Conjugation: Genie FluorRed 780
Size: 25µg, 100µg
Host Species: Mouse
Isotype: Mouse IgG2a, κ



Isotype Control: Genie Fluor Red 780 Mouse IgG2a, κ Isotype Control[C1.18.4] [Product AGEL3117]
Storage Buffer: Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Shipping: Biological ice pack at 4°C

Stability & Storage: Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

Recommended Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use. We suggest each investigator should titrate the reagent to obtain optimal results [The recommended concentration is 0.1-1 µg/10⁶ cells in 100 µL volume].
