

**Antibody Data**

<b>Product SKU:</b>	<b>AGEL3329</b>	<b>Clone:</b>	<b>B3B4</b>
<b>Applications:</b>	<b>FCM</b>		
<b>Reactivity:</b>	<b>Mouse</b>		


**Important Note:**

Centrifuge before opening to ensure complete recovery of vial contents.

**Product Information:**

**Alternate Names:** Fcεr2;Fcer2a;Lymphocyte IgE receptor;Fc-epsilon-R11 ;  
**Uniprot ID:** P20693  
**Background:** CD23 is a 45 kD protein also known as low affinity IgE Fc receptor, FcεR11, BLAST-2, Ly-42, or B6. It is a member of the Ig family, expressed on conventional B (but not B-1) cells and follicular dendritic cells. CD23 responds to high levels of IgE by downregulating IgE secretion.

**Form:** Liquid  
**Conjugation:** Genie Fluor488  
**Size:** 25µg, 100µg  
**Host Species:** Rat  
**Isotype:** Rat IgG2a, κ



**Isotype Control:** Genie Fluor 488 Rat IgG2a, κ Isotype Control[2A3] [Product AGEL3329]  
**Storage Buffer:** Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.  
**Shipping:** Biological ice pack at 4°C

**Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

**Recommended Usage:** Each lot of this antibody is quality control tested by flow cytometric analysis. Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use. We suggest each investigator should titrate the reagent to obtain optimal results [The recommended concentration is 0.1-1 µg/10<sup>6</sup> cells in 100 µL volume].