



Product Datasheet
**anti--Mouse CD134 monoclonal antibody
(FITC)**
Catalogue Code: AGIM0515

Antibody Data

| | | | |
|----------------------|----------|------------------|-------|
| Product SKU: | AGIM0515 | Clone | OX-86 |
| Applications: | FC | Conjugate | FITC |
| Reactivity: | Mouse | | |

Immunogen Information:

Host Species: Rat

Isotype: IgG1

Storage: Store at 4 °C. Avoid prolonged exposure to light.

Background: