



Product Datasheet  
**Rat CD161 NK-Cells (NKR-P1) monoclonal antibody (FITC)**  
Catalogue Code: AGIM0574

## Antibody Data

<b>Product SKU:</b>	AGIM0574	<b>Clone</b>	10-79
<b>Applications:</b>	FC	<b>Conjugate</b>	FITC
<b>Reactivity:</b>	Rat		

## Immunogen Information:

**Host Species:** Mouse

**Isotype:** IgG1

**Storage:** Store at 4 °C. Avoid prolonged exposure to light.

## Background: