



Product Datasheet
Rat CD161 NK-Cells (NKR-P1) monoclonal antibody
Catalogue Code: AGIM0575

Antibody Data

Product SKU:	AGIM0575	Clone	10-80
Applications:	FC	Conjugate	Unconjugated
Reactivity:	Rat		

Immunogen Information:

Host Species: Mouse

Isotype: IgG1

Storage: For long term storage, aliquot and freeze unused portion at -20°C in volumes appropriate for single usage. Avoid freeze/thaw cycles.

Background: