

AtH2A.Z Rabbit Polyclonal Antibody



CAB0046

Product Information

Size:

20uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

Reactivity:

Arabidopsis thaliana

Antibody Information

Recommended dilutions:

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

H2A.Z

Immunogen:

A synthetic peptide of Arabidopsis thaliana H2A.Z

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
