

CAB0077

Product Information

Product SKU:	CAB0077	Gene ID:	3833	Size:	20uL, 100uL
Clone No:	ARC1808	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	74kDa	Conjugate:	Unconjugated
Calculated MW:	74kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable microtubule binding activity and minus-end-directed microtubule motor activity. Involved in mitotic metaphase plate congression and mitotic spindle assembly. Located in membrane.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: HSET; KNSL2; KIFC1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIFC1 (Q9BW19).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.